



Stomacher® 4500 BAM



The robust Stomacher® 4500 is ideal for the BAM large volume testing procedures.

Expanding our range of Blenders focused on food laboratories requirements, the Stomacher® 4500 BAM is the perfect tool for extraction from 375g composite samples.

The Stomacher technique has become the world standard in Sample preparation. The heavy duty Stomacher 4500 BAM complements the Stomacher 400 Circulator

and provides the food laboratories with a high performance, high quality blender capable of homogenizing samples up to 4500ml.

Why the Stomacher® 4500 BAM?




Ideal for 375g samples

Efficient Large Volume
Stomacher

Volume sample size:
1000 – 4500ml

Programmable Settings

3 year Warranty

Machine Model	Paddle Type	Sample weight*	Diluent Volume*
 Stomacher 4500		375g	3375ml
Bags/Accessories	Standard Bags	Bag Rack	Bag Clip
	BA6042 BA6042/HD	BA6094	BA6097

* as per the FDA BAM protocol protocol listed 20 September 2014 :
www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods

For further information on Seward products please contact your local Seward distributor.

Seward (UK)

Europe, Middle East, Indian Subcontinent, Africa, ROW
info@seward.co.uk
+44 (0)1903 823077

Seward (Singapore)

Asia, Pacific
info@seward.co.uk
+65 6519 2283

Seward (USA)

North America, South America
info@sewardusa.com
+1 954 862 1421